

medatro® medical trolley**M01**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: $\phi 520 \times 890 \text{mm}$

Base size: $\phi 520 \times 40 \text{mm}$

Column size: $\phi 40 \times 710 \text{mm}$

Handle: $230 \times 175 \text{mm}$

Humidifier hanger: $55 \times 40 \times 16 \text{mm}$

Castor: 4pcs*4inch (with brake)

Load Capacity: 10kg

Max push speed: 2m/s

Max tilt angle: 15°

Carton dimension: $900 \times 570 \times 210 \text{mm}$

Net weight: 7kg

Gross weight: 9.6kg

