

**medatro® medical trolley****M02**

medical industry mobility solution

**Safety · Durability · Mobility · Customized LOGO**

Dimension:  $\phi$ 600\*760mm

Base size:  $\phi$ 600\*40mm

Column size:  $\phi$ 51\*610mm

Fixed bar: 400\*25\*10mm

Installing plate for fixed bar: 94\*33\*70mm

Handle: 230\*175mm

Humidifier hanger: 55\*40\*16mm

Installing plate for humidifier hanger: 55\*40\*4mm

Castor: 5pcs\*3inch (with 2brakes)

Cable hanger: 204\*65\*34mm

Load Capacity: 15kg

Max push speed: 2m/s

Max tilt angle: 15°

Carton dimension: 900\*570\*210mm

Net weight: 9.5kg

Gross weight: 12.1kg

